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Specification Sheet

siemens.com/healthcare

Excel with confidence in your hemostasis testing: Smartly designed technologies improve sample management, increase efficiency, and streamline lab workflow.

The Sysmex® CS-2500 System is a fully automated coagulation analyzer that employs smartly designed technologies such as assay-based preanalytical sample checks using PSI™ technology, automated mixing studies, automated platelet aggregation,* and clot waveform analysis† to improve sample management, increase efficiency, and streamline your lab's workflow. The system helps mid-sized to large coagulation laboratories achieve improved first-run accuracy by identifying and automatically managing potentially problematic test samples prior to analysis.

Reagent Test Menu

| | |
|-----------------------|--|
| PT | Dade® Innovin®, Thromborel® S |
| APTT | Dade Actin® Dade Actin FS, Dade Actin FSL, Pathrombin® SL |
| Fibrinogen | Dade Thrombin |
| Thrombin time | Test Thrombin, Thromboclotin* |
| Batroxobin time | Batroxobin |
| Factor deficiency | Factor II, V, VII, VIII, IX, X, XI, XII |
| Lupus anticoagulant | LA1 screening, LA2 confirmation |
| Protein C pathway | Protein C, ProC® Ac R, ProC Global, F V Leiden, INNOVANCE® Free PS Ag, Protein S Ac |
| Heparin | Berichrom® Heparin, INNOVANCE Heparin |
| Antithrombin III | Berichrom Antithrombin III (A), INNOVANCE Antithrombin |
| D-dimer | INNOVANCE D-Dimer |
| von Willebrand factor | vWF Ag, INNOVANCE vWF AC, BC von Willebrand Reagent |
| Chromogenic | Berichrom Factor VIII Chromogenic, Berichrom Factor XIII Chromogenic, Berichrom C1-Inhibitor, Berichrom α2-Antiplasmin, Berichrom Plasminogen, Berichrom Protein C |

*Sysmex CE-marked application.
†Research use only.

Sysmex CS-2500 System Specifications

Measurement

| | |
|--|---|
| Principle | Photo-optical, continuous, sequential; based on change in the transmitted light emitted from the sample with added reagent |
| Method/channels | Clotting (10 single channels), chromogenic, immunoassay, and aggregation methods |
| Measurement channels | 10 wells (mixing function using stir bars, possible for 4 wells) |
| Source lamp | Halogen lamp for measurements at 340, 405, 575, 660, and 800 nm wavelengths |
| Analysis mode | Normal and micro-sample modes |
| Analysis method | Dilution analysis, multidilution analysis (MDA), automatic reanalysis (redilution analysis, reanalysis, reflex testing), mixing testing, platelet aggregation testing, clot waveform analysis (CWA)* |
| Time resolution | Sampling can be performed for up to 1800 seconds, at intervals of 0.1 seconds |
| Measuring time | Up to 1800 seconds for each parameter |
| Incubation | 10 wells |
| Number of user-definable (open) channels | 80,000 |
| Sample Handling | |
| Type of sample | Primary tubes and/or sample cups |
| Sampling mechanism | Automated sample and standard predilution |
| Sample integrity check (PSI) | Assay-based qualitative check for hemolysis, lipemia, icterus, and proper fill of primary tube |
| Traceability of results | Operator name, test date and time, reagent lot information, reagent elapsed time onboard, HIL and sample-volume checks, test reaction position, testing and reagent table temperatures, test protocol number, dilution ratio, QC performed date, calibration curve identification, and maintenance and error logs |
| Carrier system | Continuous-access sample racks with 10-tube capacity |
| Maximum load | 50 samples, 5 racks with 10 samples per rack |
| SIA sample loading | Five priority positions |
| Storage temperature | Room temperature |
| Racks | Five sample racks, continuous loading, various specific tube holders and adapters available |
| Handling | Flexible mix of various capped and/or uncapped sample tubes and 4 mL conical sample cups |
| Primary sample probe | Liquid-level sensing, crash protection, clog detection, liquid surface verification |
| Automatic sample predilution | Samples can be diluted (0 to 19.1 to 259) for multidilution analysis (MDA) and mixing tests and can be retained for auto repeat until results are available |
| Cap piercing | Pierces caps on primary tubes |

Reagent Handling

| | | | |
|--|---|------------------------|-----------------|
| Reagent recognition | Internal reagent bar-code identification | | |
| Dispensing accuracy | Volume Dispensed | Reproducibility | Accuracy |
| Reagent probe | 50 µL | CV 3% or less | 48.0 ± 3.0 µL |
| | 100 µL | CV 1% or less | 99.0 ± 2.0 µL |
| Sample probe | 10 µL | CV 5% or less | 9.5 ± 1.0 µL |
| | 100 µL | CV 1% or less | 98.5 ± 2.0 µL |
| Volume | Reagent dispensing pipette: 20–200 µL Sample dispensing pipette: 20–250 µL (combined diluent and sample) Volumes aspirated of sample and diluent are in the range of 4–150 µL | | |
| Dispensing mechanism | Two probes: one heated for reagents and one for samples, controls, and calibrators Two catchers | | |
| Loading system | Manual placements of reagents, removable reagent trays | | |
| Reagents onboard | 40 reagent/control positions and 5 buffer/rinse solution positions | | |
| Storage temperature | 40 cooled reagent positions onboard at 10°C ± 2°C | | |
| Handling | Flexible mix of reagent positions within the reagent table, various adapters available | | |
| Mixing position | Up to 10 positions available on the reagent table | | |
| Storage capacity onboard | 1890 tests average, 3000 tests maximum | | |
| Onboard stability | System keeps reagents cooled 24 hours/day, 7 days/week, anti-evaporation caps support long onboard stability | | |
| Reagent inventory management | Tracks number of tests remaining, lot onboard stability, vial type, set date and time, and expiration date | | |
| Throughput* (tests/hour, approx.) | | | |
| PT | 180 (single-parameter analysis) | | |
| APTT | 180 (single-parameter analysis) | | |
| PT/APTT | 180 (simultaneous analysis) | | |
| PT/APTT/AT/DD | 95 (simultaneous analysis) | | |
| Operation | | | |
| Access mode | Continuous reagent-, consumable-, and sample-loading capability | | |
| Calibration | Two 12-point calibration curves with maximum of 5 repeated analyses per point and up to 10 calibration curves; one reagent lot group | | |
| Calibration curve | 250 user-defined parameters, 12 points maximum | | |
| View calibration | Graphical display of calibration curves from up to 10 different reagent lots/parameters | | |
| Auto calibration/auto QC | User-defined time interval or with new reagent vial | | |
| Temperature control | Detector 37°C ± 0.5°C, sample incubator 37°C ± 1.0°C, reagent probe 37.5°C ± 0.5°C | | |
| Sample and abnormal reaction monitoring | Coagulation curve abnormalities, overreaction to antigen, interfering assay-based HIL, sample volume check, analysis value error monitoring (checks for deviations based on mark limit values, report limits, and the results of repeated analyses) | | |

For more information, please contact your Siemens representative or visit our website.

Reaction Tubes

| | |
|------------------------|--|
| Type | Single reaction tubes |
| Loading | Automatic continuous access, 500 reaction tubes onboard |
| Reaction tube stirring | Yes |

System Fluid

| | |
|----------------------|---|
| Cleaning and rinsing | Washing solutions onboard |
| System containers | 20 L containers for water (rinse) and waste (optional) |

Computer/Printer

| | |
|-----------------------------|--|
| Workstation | PC |
| Display | 24" wide-screen LCD display with touch operation |
| Printer | Graphic printer (optional) |
| Input devices | Multitouch-compatible display (optional with monitor arm), keyboard and mouse, 2-D bar-code reader |
| Data storage | About 10,000 samples with a maximum of 60 results per sample |
| Onboard maintenance logs | Scheduled and automatically monitored routine maintenance activities via software |

Software

| | |
|----------------------|--|
| LIS interface | Sysmex CA-1000, CA-1500, and CS-2000i, CR-800, ASTM1381-95/1394-95, ASTM1381-02/1394-97 |
| Host connection | Bidirectional RS-232C serial port or via Ethernet-FCP/IP |
| Operating system | Microsoft Windows 7-based |
| System documentation | Instruction manual and reference guide online |

Power Supply

| | |
|-------------------|---|
| Operating voltage | Main unit: 100–240 V Pneumatic unit: 100–117 V/220–240 V |
| Power consumption | ≤800 VA main unit, ≤280 VA pneumatic unit (processing) |
| Main frequency | 50–60 Hz |

Environmental Conditions

| | |
|------------------------|---|
| Operating temperature | 15–30°C |
| Environmental humidity | 30–85% (no condensation except on the reagent table) |
| Atmospheric pressure | 70–106 kPa |
| Waste heat | Approx. 4000 BTU/h (1040 Kcal/h) |
| Noise level | 60 dB or below (excluding sudden noise that stops within 5 seconds and alarms) |

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Weight

| | |
|--------------------|-----------------------------|
| Main unit | Approx. 110 kg 242.5 lb |
| Pneumatic unit | Approx. 17 kg 37.5 lb |
| Terminal (PC) | Approx. 6.84 kg 15.08 lb |
| Terminal (monitor) | Approx. 5.6 kg 12.3 lb |

Quality Control

X-control, Levey-Jennings control
Multi-rule monitoring (Westgard rule)
Maximum of 1200 plots × 750 files can be saved

Compliance

| | |
|---------------------------------|---|
| Safety standards (main unit) | IEC 61010-1:2001 (2nd Edition) IEC 61010-2-081:2001 + A1:2003 IEC 61010-2-101:2002 IEC 60825-1:1993 + A1:1997 + A2:2001 EN 61010-1:2001 (2nd Edition) EN 61010-2-081:2002 + A1:2003 EN 61010-2-101:2002 EN 60825-1:1994 + A2:2001 + A1:2002 UL 61010-1:2004 (2nd Edition) CAN/CSA-C22.2 No. 61010-1:2004 (2nd Edition) CAN/CSA-C22.2 No. 61010-2-081-04 (R09) CAN/CSA-C22.2 No. 61010-2-101-04 (R09) CAN/CSA-E60825-1-03 |
|---------------------------------|---|

EMC standards

| |
|--|
| IEC 61326-1:2005 IEC 61326-2-6:2005 CISPR 11:2009/A1:2010 Group 1 Class A IEC 61000-3-2:2006/A1:2009/A2:2009 IEC 61000-3-3:2008 IEC 61326-1:2005 (Table 1) IEC 61326-2-6:2005 IEC 61000-4-2:2008 IEC 61000-4-3:2006/A1:2007/A2:2010 IEC 61000-4-4:2004/A1:2010 IEC 61000-4-5:2005 IEC 61000-4-6:2009 IEC 61000-4-8:2009 IEC 61000-4-11:2004 |
|--|

Ordering Information

| Product Description | Catalog Number |
|---------------------------|----------------|
| Sysmex CS-2500 System | 11239235 |
| Sysmex CS-2500 System kit | 11239260 |

* Sysmex CE-marked application.

* Research use only.

† Throughput values were determined by the time to first result; processing capability varies depending on the reagent used. Throughput values stated above were determined using Siemens' study protocol with PT (Thromborel S Reagent), APTT (Pathrombin SL Reagent), INNOVANCE D-Dimer Reagent, and AT (INNOVANCE AT Reagent) test applications.

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
Smartly designed technologies maximize the efficiency and productivity of your lab.

An estimated 9-15% of diagnostic errors have an impact on patient care. Many of these errors are due to the inappropriate collection, handling, or processing of samples referred for testing.

Siemens Healthcare has a proven record of providing hematology solutions that help busy labs to automate processes. In a few years, the Sysmex® CS family of systems has streamlined preanalytical sample quality checks using PS technology to help labs increase the volume of samples sent to the lab. This feature is especially important because preanalytical errors and unusable samples account for up to 70% of mistakes in the lab. In around 5% of the instances, these errors have a significant impact on patient care.

Siemens is pleased to take hematology testing to the next level with the Sysmex CS-2500 System. Built on proven PS technology employed in the Sysmex CS family of systems, the Sysmex CS-2500 System's smartly designed technologies include assay-based preanalytical sample quality checks using PS technology, automated mixing studies, automated platelet aggregation, and clot waveform analysis (CWA) for improved efficiency, operational accuracy, and reliable first-run results. Designed to increase operational efficiency and streamline workflow, the Sysmex CS-2500 System utilizes advanced technology that enables labs to:

- Raise the quality of test results
- Minimize the need for repeat testing
- Achieve cost-effective method consolidation
- Simplify operations and increase diagnostic confidence

Excel with assay-based preanalytical checks that improve sample management for accurate and consistent results.


The Sysmex CS-2500 System helps mid-volume labs achieve improved first-run accuracy by identifying and automatically managing potentially problematic test samples prior to analysis. Simultaneous multi-worklist analysis and advanced PS technologies—including assay-based hemolysis, icterus, and lipemia interference (LIL), primary tube sample volume checks and clot detection—minimize repeat testing and manual review.

- Assay-based preanalytical sample quality checks using PS technology for LIL.

Interference provides accurate and reliable results on the first run.

- Simultaneous multi-worklist scanning of clotting reactions at 340, 405, 575, 640, and 800 nm helps to reduce the effects of interfering substances by automatically selecting optimal wavelengths.
- Ten flexible reaction detectors enable high-capacity performance for a variety of test profiles to maximize productivity of the lab.


The level of hemolysis causing clinically important changes is test-specific; therefore, test-specific thresholds should be used to assess patient samples that have hemolyzed during collection and/or processing.

Sysmex CS is the application. Please read or refer to the effect of hemolysis of hematology test. Refer to the attached LIL, icterus and lipemia interference guide for more information.

1 ks – Koagulometer: analizátor SIEMENS - Sysmex CA-1500

The Sysmex® CA-1500 System quickly analyzes large quantities of samples with a high-level of accuracy. The compact, fully automated system performs clotting, chromogenic, and immunologic tests. Its advanced technology allows accurate sample identification and multiple dilution analysis (MDA) for detecting inhibitors during factor assays.

- Optional second-generation cap piercing enhances safety and productivity
- Accurate sample identification with bar-code reading at the time of aspiration
- Comprehensive menu of clotting, chromogenic, and immunologic tests, including routine tests and the most commonly ordered specialty tests
- A rapid turnaround environment—automated D-dimer results within 7 minutes
- Simultaneous calibration and routine analysis
- Bar-code reader accurately identifies STAT patient samples and reagents

The Sysmex® CA-1500 System, with optional cap-piercing, is ideal for mid-volume laboratories requiring an easy to use analyzer capable of running clotting, chromogenic, and immunologic assays simultaneously. The system features multiple dilution analysis (MDA) for detecting inhibitors during factor assay analysis and accurate sample identification. By offering true 15-parameter random access testing and walk-away capability, the system offers both efficiency and flexibility for your laboratory.

Increases productivity with rapid turnaround time

- Approximately 120 PT tests per hour
- Approximately 80 PT/APTT simultaneous tests per hour
- Detects inhibitors during factor assay testing with multiple dilution analysis (MDA)

Streamlines workflow through intelligent calibration and loading features

- Saves technicians' valuable time through simultaneous calibration and routine analysis
- Features the ability to load multiple bottles or multiple lots of reagents, decreasing preparation time
- Increases technicians' productivity with user-defined auto-repeat, re-dilution, and reflex testing

Provides confidence in patient results with safe and accurate sample-handling features

- Enhances safety and productivity with second-generation cap-piercing (optional)
- Eliminates manual input and increases safety with bar-code reader for positive STAT ID of patient samples and reagents
- Ensures correct sample results by accurately identifying samples with bar-code reading at the time of aspiration

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