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Error: Mst – 08 – XXX* – XXX* PT12 Short-circuit <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>Short-circuited wire of the temperature sensor</td><td>Service</td><td>Check, replace the sensor</td></tr> <tr> <td>Inaccurate setting of pressure and temperature sensors</td><td>Service</td><td>Check, calibration</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	Short-circuited wire of the temperature sensor	Service	Check, replace the sensor	Inaccurate setting of pressure and temperature sensors	Service	Check, calibration			
Cause	Removal by	Action to perform												
Short-circuited wire of the temperature sensor	Service	Check, replace the sensor												
Inaccurate setting of pressure and temperature sensors	Service	Check, calibration												
Error: Mst – 09 – XXX* – XXX* PE12 Disconnected <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>Disconnected wire of the temperature sensor</td><td>Service</td><td>Check, replace the sensor</td></tr> <tr> <td>Inaccurate setting of pressure and temperature sensors</td><td>Service</td><td>Check, calibration</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	Disconnected wire of the temperature sensor	Service	Check, replace the sensor	Inaccurate setting of pressure and temperature sensors	Service	Check, calibration			
Cause	Removal by	Action to perform												
Disconnected wire of the temperature sensor	Service	Check, replace the sensor												
Inaccurate setting of pressure and temperature sensors	Service	Check, calibration												
Error: Mst – 10 – XXX* – XXX* PE12 Short-circuit <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>Short-circuited wire of the temperature sensor</td><td>Service</td><td>Check, replace the sensor</td></tr> <tr> <td>Inaccurate setting of pressure and temperature sensors</td><td>Service</td><td>Check, calibration</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	Short-circuited wire of the temperature sensor	Service	Check, replace the sensor	Inaccurate setting of pressure and temperature sensors	Service	Check, calibration			
Cause	Removal by	Action to perform												
Short-circuited wire of the temperature sensor	Service	Check, replace the sensor												
Inaccurate setting of pressure and temperature sensors	Service	Check, calibration												
Error: Mst – 12 – XXX* – XXX* Slv – no response <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>Automatics error</td><td>Operator</td><td>Wait for 10 s, start the cycle</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	Automatics error	Operator	Wait for 10 s, start the cycle						
Cause	Removal by	Action to perform												
Automatics error	Operator	Wait for 10 s, start the cycle												
Error: Slv – 02 – XXX* – XXX* Door unblocked <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>The safety device registered improper closing of the door</td><td>Service</td><td>Adjust the door switches.</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	The safety device registered improper closing of the door	Service	Adjust the door switches.						
Cause	Removal by	Action to perform												
The safety device registered improper closing of the door	Service	Adjust the door switches.												
Error: Slv – 06 – XXX* – XXX* Generator interrupted <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>Pump fault</td><td>Service</td><td>Check the pump</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	Pump fault	Service	Check the pump						
Cause	Removal by	Action to perform												
Pump fault	Service	Check the pump												
Error: Slv – 07 – XXX* – XXX* Vacuum pump blocked <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>El. vacuum pump protection registered the current rise above the max. limit</td><td>Service</td><td>Check, calibration</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	El. vacuum pump protection registered the current rise above the max. limit	Service	Check, calibration						
Cause	Removal by	Action to perform												
El. vacuum pump protection registered the current rise above the max. limit	Service	Check, calibration												
Error: Slv – 08 – XXX* – XXX* Low temperature <table border="1"> <thead> <tr> <th>Cause</th><th>Removal by</th><th>Action to perform</th></tr> </thead> <tbody> <tr> <td>Leakage</td><td>Operator</td><td>Perform the vacuum test Restart the cycle</td></tr> <tr> <td>Incorrect setting of sterilization cycle parameters.</td><td>Service</td><td>Check the parameters</td></tr> <tr> <td>Inaccurate setting of pressure and temperature sensors</td><td>Service</td><td>Check, calibration</td></tr> </tbody> </table>			Cause	Removal by	Action to perform	Leakage	Operator	Perform the vacuum test Restart the cycle	Incorrect setting of sterilization cycle parameters.	Service	Check the parameters	Inaccurate setting of pressure and temperature sensors	Service	Check, calibration
Cause	Removal by	Action to perform												
Leakage	Operator	Perform the vacuum test Restart the cycle												
Incorrect setting of sterilization cycle parameters.	Service	Check the parameters												
Inaccurate setting of pressure and temperature sensors	Service	Check, calibration												

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
* Does not occur in this program.

P7 – Customer program, whose parameters are specified according to customer's requirements on program course.

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